

# **Anti-Oxytocin Receptor Rabbit Monoclonal Antibody**

Catalog # ABO15661

## **Specification**

#### Anti-Oxytocin Receptor Rabbit Monoclonal Antibody - Product Information

Application WB, IHC
Primary Accession P30559
Host Rabbit
Isotype IgG

Reactivity Rat, Human, Mouse

Clonality Monoclonal Format Liquid

**Description** 

Anti-Oxytocin Receptor Rabbit Monoclonal Antibody . Tested in WB, IHC application. This antibody reacts with Human, Mouse, Rat.

## Anti-Oxytocin Receptor Rabbit Monoclonal Antibody - Additional Information

Gene ID 5021

**Other Names** 

Oxytocin receptor, OT-R, OXTR

Calculated MW 43 kDa KDa

**Application Details** 

WB 1:500-1:2000<br>IHC 1:50-1:200

#### Contents

Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol, 0.4-0.5mg/ml BSA.

#### **Immunogen**

A synthesized peptide derived from human Oxytocin Receptor

**Purification** 

Affinity-chromatography

Storage Store at -20°C for one year. For short term

storage and frequent use, store at 4°C for

up to one month. Avoid repeated

freeze-thaw cycles.

## **Anti-Oxytocin Receptor Rabbit Monoclonal Antibody - Protein Information**

**Name OXTR** 



#### **Function**

Receptor for oxytocin. The activity of this receptor is mediated by G proteins which activate a phosphatidylinositol-calcium second messenger system.

#### **Cellular Location**

Cell membrane; Multi-pass membrane protein.

### **Anti-Oxytocin Receptor Rabbit Monoclonal Antibody - Protocols**

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

### **Anti-Oxytocin Receptor Rabbit Monoclonal Antibody - Images**




Figure 1. Western blot analysis of OXTR using anti-OXTR antibody (M01566).

Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions.

Lane 1: human Jurkat whole cell lysates,

Lane 2: human PC-3 whole cell lysates,

Lane 3: human A549 whole cell lysates,

Lane 4: human Caco-2 whole cell lysates,

Lane 5: rat PC-12 whole cell lysates,

Lane 6: mouse RAW264.7 whole cell lysates.

After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-OXTR antigen affinity purified monoclonal antibody (Catalog # M01566) at 1:500 overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:500 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog #



EK1002) with Tanon 5200 system. A specific band was detected for OXTR at approximately 43 kDa. The expected band size for OXTR is at 43 kDa.




Figure 2. IHC analysis of OXTR using anti-OXTR antibody (M01566).

OXTR was detected in a paraffin-embedded section of rat ovary tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 1:50 rabbit anti-OXTR Antibody (M01566) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.




Figure 3. IHC analysis of OXTR using anti-OXTR antibody (M01566).

OXTR was detected in a paraffin-embedded section of rat ovary tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 1:50 rabbit anti-OXTR Antibody (M01566) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.