

Anti-CD63 (Extracellular region) Antibody
Catalog # AN1701**Specification****Anti-CD63 (Extracellular region) Antibody - Product Information**

Application	WB
Primary Accession	P08962
Host	Mouse
Clonality	Mouse Monoclonal
Isotype	IgG2a
Calculated MW	25637

Anti-CD63 (Extracellular region) Antibody - Additional Information

Gene ID 967

Other Names

CD63 antigen, Granulophysin, Lysosomal-associated membrane protein 3, LAMP-3, Melanoma-associated antigen ME491 OMA81H, Ocular melanoma-associated antigen, Tetraspanin-30, Tspan-30, MLA1, TSPAN30

Target/Specificity

Tetraspanins comprise a large superfamily of cell surface-associated membrane proteins with four transmembrane domains. On cell membranes, tetraspanins form networks of various proteins called tetraspanin-enriched microdomains (TEMs). CD63 was the first characterized tetraspanin and it is found in TEMs, as well as late endosomes and lysosomes. In late endosomes, CD63 is enriched on the intraluminal vesicles, and can be secreted as exosomes through fusion of endosomes with the plasma membrane. The complex localization pattern of CD63 suggests that its intracellular trafficking and distribution must be tightly regulated. CD63 contains N-link glycosylation sites that produce diverse CD63 molecules that range from 30 to 60 kDa depending on cell type. CD63 is an important exosomal marker in cancer cells, and may be involved in cancer progression and metastasis.

Format

Protein G Purified

Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Precautions

Anti-CD63 (Extracellular region) Antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Shipping


Blue Ice

Anti-CD63 (Extracellular region) Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

Anti-CD63 (Extracellular region) Antibody - Images



Anti-CD63 (Extracellular region) Antibody - Background

Tetraspanins comprise a large superfamily of cell surface-associated membrane proteins with four transmembrane domains. On cell membranes, tetraspanins form networks of various proteins called tetraspanin-enriched microdomains (TEMs). CD63 was the first characterized tetraspanin and it is found in TEMs, as well as late endosomes and lysosomes. In late endosomes, CD63 is enriched on the intraluminal vesicles, and can be secreted as exosomes through fusion of endosomes with the plasma membrane. The complex localization pattern of CD63 suggests that its intracellular trafficking and distribution must be tightly regulated. CD63 contains N-link glycosylation sites that produce diverse CD63 molecules that range from 30 to 60 kDa depending on cell type. CD63 is an important exosomal marker in cancer cells, and may be involved in cancer progression and metastasis.