

**MOUSE IgG3 isotype control Allophycocyanin**  
**Monoclonal IgG3 , Allophycocyanine**  
**Catalog # ASR3324****Specification**

---

**MOUSE IgG3 isotype control Allophycocyanin - Product Information**

Description	<b>MOUSE IgG3 isotype control Allophycocyanin Conjugated</b>
Conjugate	<b>Allophycocyanine</b>
FP Value	<b>1-2 moles Allophycocyanine per mole of Mouse IgG3</b>
Clonality	<b>Monoclonal</b>
Application	<b>,4,10,</b>
Application Note	<b>ELISA 1:2000-1:20,000;FlowCytometry 1:1000-1:5000</b>
Physical State	<b>Liquid (sterile filtered)</b>
Host Isotype	<b>IgG3</b>
Species of Origin	<b>Mouse</b>
Stabilizer	<b>None</b>
Preservative	<b>0.01% (w/v) Sodium Azide</b>

**MOUSE IgG3 isotype control Allophycocyanin - Additional Information****Shipping Condition**

Wet Ice

**Purity**

Mouse Isotype control has been prepared from concentrated cell culture supernatant by immunoaffinity chromatography using protein A. In an Ouchterlony double diffusion assay the material is non-reactive with antisera to mouse IgG1, IgG2a, IgG2b, IgM, and IgA. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Mouse IgG and anti-Mouse serum. Light and heavy chain composition has been confirmed.

**Storage Condition**

Store vial at 4° C prior to opening. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. DO NOT FREEZE. This product is light sensitive.

**Precautions Note**

This product is for research use only and is not intended for therapeutic or diagnostic applications.

**MOUSE IgG3 isotype control Allophycocyanin - Protein Information****MOUSE IgG3 isotype control Allophycocyanin - Protocols**

Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

#### **MOUSE IgG3 isotype control Allophycocyanin - Images**

#### **MOUSE IgG3 isotype control Allophycocyanin - Background**

Mouse isotype controls are used in flow cytometry, western blot and ELISA and differentiate between immunoglobulin classes and subclasses. Isotype controls allow for the genetic variations or differences in the constant regions of the heavy and light chains. In mouse there are six relevant heavy chain isotypes and two light chain isotypes: heavy chain  $\alpha$  - IgA,  $\gamma$  - IgG 1, 2a, 2b, 3 and  $\mu$  - IgM, light chain  $\kappa$  and  $\lambda$ .